Strain Seer - 2D Strain Analysis Tool

■ What is 2D Strain Analysis?

2D strain analysis measures material deformation using 5 fiducial markers and 2 scale markers to calculate:

- εxx: Stretching/compression in x-direction
- εyy: Stretching/compression in y-direction
- Exy: Shearing deformation

How It Works


- 1. Mark 5 fiducial points (red) and 2 scale points (green)
- 2. Enter scale length to convert pixels to real-world units
- 3. Input deformation distance for each image
- 4. View strain analysis and regression results

Point Labeling Convention

```
0 ----- 1
| \ / |
| 4 |
| / \ |
3 ----- 2
```

- Points 0-3: Corner points (clockwise from top-left)
- Point 4: Center point
- Scale points: Calibration markers (green)

№ Why We Need Scale Points



Scale points ensure:

- Consistent measurements across images
- Accurate pixel-to-real-world conversion
- Reproducible results

🖋 Step-by-Step Guide

- 1. Upload Images
 - Supported: PNG, JPG, JPEG, GIF, BMP
 - o 🛕 Important: Upload files in order
 - Order determines sequence in analysis and plots
- 2. Annotate Points
 - 5 red fiducial points
 - 2 green scale points
 - Enter deformation distance

https://strain-seer.streamlit.app

3. Run Analysis • Set scale length • View results Export data **Tips for Best Results** Clear fiducial markers Known scale distance Well-lit, focused images Correct point order (0-4) Upload images in sequence Choose image files Drag and drop files here Browse files Limit 200MB per file • PNG, JPG, JPEG, GIF, BMP deformation-30mm-force.png 1.5MB deformation-25mm-force.png 1.6MB deformation-20mm-force.png 1.3MB × Files Status Fiducial Points 1 deformation-00mm-force.png 0.000000 5.000000 2 deformation-05mm-force.png 3 deformation-10mm-force.png 10.000000 15.000000 4 deformation-15mm-force.png 20.000000 5 deformation-20mm-force.png 6 deformation-25mm-force.png 25.000000 7 deformation-30mm-force.png 30.000000 Select image to annotate deformation-30mm-force.png

https://strain-seer.streamlit.app



Point Coordinates

Туре	Order	Х	
	1	299.5	495.5
	2	535.4	522.7
	3	508.2	733.1
	4	201.1	729.7
	5	381.0	656.7
	6	703.4	52.6
	7	688.1	687.3

Deformation Distance

Enter deformation distance (mm)

30.00

- +

Strain Analysis



Next steps:

- 1. Set scale length
- 2. View results
- 3. Export data


♦ Scale Calibration


Enter physical length of scale markers (mm)

Scale Length (mm)


50.00


∠ X-Axis Strain (εxx)






https://strain-seer.streamlit.app 3/4







Regression Equation:

 $\epsilon xx = 1.27e-17 \text{ mm}^{-1}x + 2.47e-32$

R² Value: 0.121212 P-value: 0.444113

Distance (mm)	εхх
0.00	-2.22e-16
5.00	0.00e+00
10.00	2.22e-16
15.00	4.44e-16
20.00	2.22e-16
25.00	8.88e-16
30.00	-2.22e-16

Regression Equation:

εyy = 1.59e-17 mm⁻¹x + -2.06e-16

R² Value: 0.211864 P-value: 0.298667

Distance (mm)	εγγ
0.00	-4.44e-16
5.00	-2.22e-16
10.00	0.00e+00
15.00	2.22e-16
20.00	6.66e-16
25.00	2.22e-16
30.00	-2.22e-16

Regression Equation:

 $\varepsilon xy = -5.05e-18 \text{ mm}^{-1}x + 5.20e-17$

R² Value: 0.389734 **P-value:** 0.133993

Distance (mm)	εχγ
0.00	2.02e-18
5.00	3.50e-17
10.00	2.83e-17
15.00	3.20e-17
20.00	8.94e-18
25.00	-2.05e-16
30.00	-6.70e-17

Data Export

Raw Data

Available formats:

- JSON: File data + scale length
- CSV: Strain components

Download Raw Data (JSON)

Download Raw Data (CSV)

Analysis Results

Available formats:

- JSON: Regression parameters
- CSV: Regression data

Download Analysis Results (JSON)

Download Analysis Results (CSV)

⚠ New Analysis

To start fresh:

- 1. Refresh page
- 2. Upload images
- 3. Mark points
- 4. Set scale length